

[PAR Scientific A/S](#)[Products](#)[Support](#)[International contacts](#)[What's new](#)[.: Home -> Products -> Fixation -> VacFix Accessories:.](#)

VacFix[®] Accessories



Vacuum Pumps.

Par Scientific offers two handy vacuum pumps for evacuating the cushions. Both models are cordless providing easy operation. Selection of type depends on demands.



VacFix[®] Storage Carts

The cart is available in 2 standard sizes to meet the most basic demands.

VacFix[®] Storage Cart increases the cushion life and decreases the demand for additional space in each accelerator room.

The carts are equipped with wheels for easy transportation.



Knee Pillows

Knee pillows are an essential part of a good Pelvic/Prostate fixation.

Providing an excellent patient comfort and ensures a good contact to the vacuum cushion and a solid fixation.



Fixation Table

The VacFix[®] Fixation Table is designed to provide adjustment flexibility and stability yet being a light weight system that is easy to move. The large size table top is constructed in co-operation with skilled VacFix[®] users world wide to meet any clinical demand during mould.



Marking Tape

Marking tape in two different colours are used in a variety of clinical situations as a vital helping tool. Being adhesive on one side it can be used anywhere where accurate removable markings are required.

Please investigate detailed description and use in the clinical fixation studies.



Position Guides

Position Guides (male and female) are used on simulator table top and accelerator table top in identical patterns allowing the associated guides positioned at the back of each cushion to match exactly thereby obtaining exact and identical position of the cushion.

Velcro Tape



Velcro tape - or hook and velour tape is used mainly to mount thermoplastic sheet to the vacuum cushion. Please investigate detailed description and use in the clinical fixation studies.

© PAR Scientific A/S. Optimized for 800x600. Updated: Friday January 3 2003 20⁵⁰. | [E-mail](#) | [Top](#) | [Go back](#)